

HOME > TERMS >
 ABOUT > HELPS/FAQS >
 SPONSORS > NEWSLETTER >

MATERIALS > APPLICATIONS > INDUSTRIES >
 CLASSIFIEDS > COURSES > EVENTS >
 JOBS > NANO JOURNAL > NEWS >

ADVERTISE ON AZONANO

ADVERTISE HERE


EDS & WDS X-ray microanalysis systems

Detailed structural and chemical information at the micro and nanoscale



Cenamps

Applying small-scale technologies

Nanotechnology News

[+ Add to Google](#)

[Email to a Friend](#)

[back](#)



PI

Nanopositioning & Scanning Systems

New Single-Molecule Imaging System

Scientists are able to view active molecules within a biological motor of the nanometer scale with the help of a new imaging system far more sensitive and powerful than existing optical microscopes.

Raman Spectrometers



HORIBA JOBIN YVON

A [Purdue University](#) researcher has created a single-molecule imaging system to view deoxyribonucleic acid (DNA), ribonucleic acid (RNA) and other tiny biological molecules. This ability helps settle a seven-year debate within the virology and nanomedicine fields over the shape and structure of a tiny biological motor that has potential use in nanotechnology and nanomedicine, including the diagnosis and treatment of diseases such as cancer, AIDS and influenza.

Special Issue



Applied Physics A
Materials Science & Processing

Semiconductor Nanowires.

Guest Ed: J. E. Fischer

Scientists had disputed the number of packaging ribonucleic acid (pRNA) molecules contained in the DNA-packaging motor of the phi29 virus. The number of these molecules present determines the shape of the motor and expands understanding of the way it works. The new imaging system definitively concludes that six pRNA molecules were present. The research, funded by the National Institutes of Health, will be published in the upcoming issue of the European Molecular Biology Organization Journal, EMBOJ.

"The microscopes commercially available were not sensitive enough to distinguish individual molecules of this scale," said

Ads by Google

[Proton Beam Therapy](#)
 See how Proton Beam Therapy can treat cancer. Watch the true story.
www.hitachi.com/truestories

[Nikon Confocal](#)
 New Family of Confocal Microscope Systems from Nikon
www.nikonconfocal.com

[Nanotech Research Pioneer](#)
 Learn More About Ecology Coatings' 100% Solids UV Curable Technology
www.ecologycoatings.com

[Medical Imaging Research](#)
 High Quality Diagnostic Imaging & Patient Comfort From Siemens.
www.usa.siemens.com/Medical

[Hybrid LCMS System](#)
 Hitachi NanoFrontier LCMS Ultimate in Hybrid Technology
www.hitachi-hta.com

[Next Generation AFM/SPMs](#)
 Advanced atomic force microscopes for nanotechnology applications
www.AsylumResearch.com

[Imaging](#)
 Auto Image Acquisition & Analysis Improve Output w/Speed & Efficiency
www.moleculardevices.com

[NanoAndMore USA](#)
 AFM tips from budget to research imaging software - special offers
www.nanoandmore.com

ADVERTISE HERE

Peixuan Guo, professor of molecular virology and biomedical engineering and creator of the imaging system. "Our system, through its highly sensitive detection system and dual color viewing ability, is capable of distinguishing and counting individual molecules within a nanodevice."

Nanotechnology deals with materials within the nanometer scale. A nanometer is one-billionth of a meter, or 100,000 times smaller than the diameter of a human hair. The most powerful optical microscopes available offer a resolution of 200 nanometers, meaning anything closer together than 200 nanometers would blur together and appear as a single unit.

"Many of the molecules involved in the cellular processes of the human body are much smaller than this 200-nanometer diffraction limit," Guo said. "For example, the phi29 motor is approximately 20 nanometers in size and a strand of DNA is only 2 nanometers wide."

Electron microscopes, which use a beam of electrons instead of light to detect a specimen, can reach a resolution of tenths of a nanometer. However, materials that have a versatile or flexible structure but are not electron-dense, such as RNA, are not easily imaged using electron microscopes, he said.

Active molecules performing their functions can be observed in real time through Guo's system, which can detect the motion of the molecules and also count individual molecules within a few tens of nanometers.

"This was greatly needed for nanotechnology research," he said. "Researchers need a clear view of their subjects to fully understand them. Direct observation of these molecules will answer many of the questions facing this field and will help in the design of novel nanodevices. The strength of our system is the sensitivity of detection."

Guo, who leads Purdue's NIH-funded National Nanomedicine Development Center and also is a member of Purdue's Cancer Center, used readily available parts to construct the innovative imaging system.

The system, called a single-molecule dual viewing total internal reflection fluorescence imaging system (SMDV-TIRF) is an improvement of total internal reflection microscopy. Its improved sensitivity comes from the use of a newly assembled laser combiner to control the release of two or more lasers as the light source and the direct delivery of this laser light through optical fibers to the detector.

Other modifications, such as an advanced cooling system and the placement of the quartz prism on the top of the chamber, minimize background noise and leakage of light that could confuse the detector. In addition, the application of binomial distribution for dual view technology has led to new findings in imaging.

The molecules of interest are marked with red and green fluorescent labels that emit signals. These signals are divided into different paths according to their colors before they reach the detector and a computer analyzes the information and converts the signals to images.

The modifications to the system assure that the signal from the labels retain their strength and are clearly detected, Guo said.

The system recognizes when different labels overlap, a sign of interaction

between the different types of molecules, and presents a yellow color to differentiate it from the red and green markers. As the system analyzes the sample, the computer displays an image and a graph.

The image displays dots of red, green and yellow, showing the location of the DNA and RNA and where and when the two overlap.

The graph displays lines similar to stairsteps, which are used to count the individual molecules. As the light shines on the fluorescent markers to illuminate them, it also degrades and, eventually, destroys the marker. As this destruction, called photobleaching, occurs, less light is sent to the detector.

"The fluorescent molecules naturally depart one by one when exposed to intensive illumination," Guo said. "The graph drops as each molecule's marker fades until none are left and the graph reaches the baseline. Each stairstep on the final graph represents one single molecule with one marker."

Guo and his team, including co-authors Hui Zhang and Dan Shu, and postdoctoral researcher Wulf-Dieter Moll were able to count the number of pRNA molecules present and observed the motion of the phi29 motor as it packaged DNA. By observing the interaction of the DNA and RNA they were able to determine the structure and mechanics of the nanomotor.

"The pRNA molecules are unstable and quickly degrade," Guo said. "This caused confusion over the number of pRNA molecules present. Researchers were unable to directly count the molecules when the motor is in motion and had to estimate the number of molecules present based on the chemical composition of large samples. The degradation of the RNA affected the estimates, and no one could solve the problem."

Guo said to construct this novel system, he received assistance and advice of many single molecule biophysicists, including Nils Walter, David Rueda, Peter Stockley, Taekjip Ha, Eckhard Jankowsky, Toshio Yanagida, Chris Meiners, Meredith Lambert, Faqing Huang and Mark Browne from Ando Technology.

Guo's team next plans to count the individual protein components within the phi29 DNA packaging motor and will further characterize the real-time motion of the motor.

<http://www.purdue.edu>

Posted 31st January 2007

[Nanotechnology News Archive](#)

[↑](#)
top
[↩](#) back



AZoNano is grateful for the support provided by our [sponsors](#) to both AZoNano.com and to the authors and peer reviewers of [AZoJono - Journal of Nanotechnology Online](#) - open access to leading Nanotech Science.

Please visit our [sponsors page](#) or click the links below to view the world leading products and services our [sponsors](#) have to offer. To discuss an AZoNano sponsor package please contact our [Sales and Marketing Manager](#)



Other AZoNetwork Sites | AZoM.com | AZoBuild.com | News-Medical.Net

[HOME >](#)
[ABOUT >](#)
[SPONSORS >](#)

[TERMS >](#)
[HELPS/FAQS >](#)
[NEWSLETTER >](#)

MATERIALS >	APPLICATIONS >	INDUSTRIES >
CLASSIFIEDS >	COURSES >	EVENTS >
JOBS >	NANO JOURNAL >	NEWS >

AZoNano - The A to Z of Nanotechnology...Copyright © 2007 by AZoM.com Pty.Ltd